


## WP 300.05

### Compression plates for compression tests, large



#### Learning objectives/experiments

- pressure experiments on the WP 300 experimental unit for own specimens

#### Specification

- [1] accessory set of pressure plates for own specimens on the WP 300 experimental unit

#### Technical data

- 2 pressure plates
  - LxW: 160x60x50mm
  - material: steel, hardened

LxWxH: 2x 160x60x50mm  
Weight: approx. 8kg

#### Scope of delivery

- 2 pressure plates
- 1 set of accessories

#### Description

- **easy to install sample holder for pressure testing of own specimens on the WP 300 experimental unit**

This accessory for WP 300 makes it possible to carry out pressure experiments. The compressive strength is determined as an important characteristic of a material. Furthermore, the compression yield point can be determined as the onset of plastic deformation of a material, at which no compressive strength can be determined.

The accessory contains two pressure plates between which a specimen is positioned. The assembled compression mechanism is installed in the pressure section of the WP 300 experimental unit, between the lower crossmember and the crosshead.

In experiments, a uniaxial state of stress is produced in a geometrically defined specimen. This state of stress is produced by an external load on the specimen in the longitudinal direction via a compressive force. Then a uniform normal stress distribution prevails in the test cross-section of the specimen.

In order to determine the strength of the material, the load on the specimen is slowly and steadily increased, until the specimen ruptures.

Materials with a relatively low compressive strength or different geometry that require a larger contact surface, can be studied e.g. assembly foam, cardboard boxes, plastic bottles (lab-own compression specimens).

## **WP 300.05**

### **Compression plates for compression tests, large**

Required accessories

WP 300

Materials testing, 20kN