

HM 240 Principles of air flow



Description

wide range of accessories for basic experiments with air flow
record a fan characteristic

HM 240 is part of a series that allows experiments on the fundamentals of air flow. The software for data acquisition and visualisation makes the experiments especially clear and enables fast execution of experiments with reliable results.

The experimental unit includes a radial fan, which can be used to generate flow velocities up to 9m/s. An inlet contour on the intake side ensures a low-turbulence flow and thus a homogeneous velocity distribution in the measuring section. A throttle valve on the end of the pressure pipe can be used to adjust the air flow to allow the fan characteristic curve to be recorded. When used in conjunction with the power meter HM 240.02 it is possible to determine the efficiency of the fan.

Further accessories for additional experiments can be attached in the intake pipe: electronic total pressure sensor HM 240.03, pressure distribution on a cylinder HM 240.04 and heat transfer at a cylinder in transverse flow HM 240.06. To study the friction losses, the intake pipe is replaced with pipe elements from HM 240.05 (straight pipes, pipe bends and pipe angles).


Measuring points along the measuring section allow temperature, pressure and velocity measurements to be taken. The flow rate is determined by means of the inlet contour and the pressure measurement. The measured values are transmitted directly to a PC via USB. The data acquisition software is included.

Learning objectives/experiments


- recording a fan characteristic
- in conjunction with the power meter HM 240.02
 - determining the fan efficiency
- in conjunction with corresponding accessories
 - velocity distribution in the pipe
 - velocity distribution behind a cylinder subject to transverse incident flow
 - pressure distribution around a cylinder subject to transverse incident flow
 - friction losses in pipes, pipe bends and pipe angles
 - recording the cooling curve of a copper cylinder subject to incident flow
 - determining the heat transfer coefficients from the cooling curve



HM 240 Principles of air flow




1 air outlet, 2 throttle valve for adjusting the air flow, 3 measuring point for temperature, 4 fan, 5 switch box with pressure transducer, 6 delivery pipe, 7 measuring point for pressure, 8 connector for accessory HM 240.04 / HM 240.06, 9 air inlet, 10 intake pipe, 11 connector for Pitot tube HM 240.03



Representation of a fan characteristic

blue: measured values, red fan characteristic; p pressure, Q volumetric flow rate



Software screenshot

Specification

- [1] investigation of the principles of air flow
- [2] transparent intake pipe with mounting options for additional accessories
- [3] inlet contour minimises turbulence on the intake side
- [4] throttle valve on the delivery pipe to adjust the air flow
- [5] electronic measurement of temperature and pressure
- [6] determine velocity by means of the dynamic pressure
- [7] determine flow rate via differential pressure
- [8] GUNT software for data acquisition via USB under Windows 10

Technical data

Radial fan

- max. power consumption: 90W
- speed: 2800min⁻¹-
- max. flow rate: 460m³/h
- max. differential pressure: 480Pa

Delivery pipe

- outer Ø: 110mm
- inner Ø: 99,4mm

Intake pipe

- outer Ø: 140mm
- ∎ inner Ø: 134,4mm

Measuring ranges

- pressure: 1x ±10mbar
- pressure: 2x ±1mbar
- temperature: 0...200°C

230V, 50Hz, 1 phase 230V, 60Hz, 1 phase; 120V, 60Hz, 1 phase UL/CSA optional LxWxH: 850x450x600mm Weight: approx. 23kg

Required for operation

PC with Windows

Scope of delivery

- 1 experimental unit
- 1 set of hoses
- 1 GUNT software + USB cable
- 1 set of instructional material

G.U.N.T. Gerätebau GmbH, Hanskampring 15-17, D-22885 Barsbüttel, Telefon (040) 67 08 54-0, Fax (040) 67 08 54-42, Email sales@gunt.de, Web www.gunt.de We reserve the right to modify our products without any notifications. Page 2/3 - 11.2023



HM 240 Principles of air flow

Optional accessories

for Remote Learning]
GU 100	Web Access Box
with	
HM 240W	Web Access Software

Other accessoriesHM 240.02Power meterHM 240.03Electronic total pressure sensorHM 240.04Pressure distribution on a cylinderHM 240.05Pressure losses in pipe elementsHM 240.06Heat transfer at a cylinder in transverse flowWP 300.09Laboratory trolley

G.U.N.T. Gerätebau GmbH, Hanskampring 15-17, D-22885 Barsbüttel, Telefon (040) 67 08 54-0, Fax (040) 67 08 54-42, Email sales@gunt.de, Web www.gunt.de We reserve the right to modify our products without any notifications. Page 3/3 - 11.2023